THC Rapid Test Cassette (Oral fluid)

Package Insert

REF DTH-P802
(382310)

English

A rapid test for the qualitative detection of Marijuana in human Oral Fluid.
For professional *in vitro* diagnostic use only.

[INTENDED USE]

The THC Rapid Test Cassette (Oral fluid) is a rapid chromatographic immunoassay for the detection of Marijuana in human oral fluid at a cut-off concentration of 3.5 ng/mL. This assay provides only a preliminary test result. A more specific alternate chemical method must be used in order to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

[SUMMARY]

11-nor- Δ^9 -tetrahydrocannabinol-9-carboxylic acid (Δ^9 -THC-COOH), the metabolite of THC (9-tetrahydrocannabinol), is detectable in oral fluid shortly after use. The detection of the drug is thought to be primarily due to the direct exposure of the drug to the mouth (oral and smoking administrations) and the subsequent sequestering of the drug in the buccal cavity. Historical studies have shown a window of detection for THC in oral fluid of up to 14 hours after drug use.

[PRINCIPLE]

The THC Rapid Test Cassette (Oral fluid) is a rapid chromatographic immunoassay based on the principle of competitive binding. Drugs that may be present in the oral fluid specimen compete against their respective drug conjugate for binding sites on their specific antibody. During testing, a portion of the oral fluid specimen migrates upward by capillary action. Marijuana, if present in the oral fluid specimen below 3.5 ng/mL, will not saturate the binding sites of the antibody coated particles in the cassette. The antibody coated particles will then be captured by immobilized THC conjugate and a visible colored line will show up in the test line region. The colored line will not form in the test line region if the Marijuana level is above 3.5 ng/mL because it will saturate all the binding sites of anti-Marijuana antibodies. A drug-positive oral fluid specimen will not generate a colored line in the test line region because of drug competition, while a drug-negative oral fluid specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

[REAGENTS]

The test contains rabbit monoclonal THC antibody-coupled particles and THC-protein conjugate. A goat polyclonal antibody against gold-protein conjugate at the control line, and a dye pad which contains colloidal gold particles coated with rabbit monoclonal antibody specific to Marijuana.

[PRECAUTIONS]

- Do not use after the expiration date.
- The test should remain in the sealed pouch until use.
- Oral fluid is not classified as biological hazard unless derived from a dental procedure.
- The used collector and cassette should be discarded according to local regulations.

[STORAGE AND STABILITY]

Store as packaged in the sealed pouch at 2-30°C. The test is stable through the expiration date printed on the sealed pouch. The test cassettes must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

[SPECIMEN COLLECTION AND PREPARATION]

The oral fluid specimen should be collected using the collector provided with the kit. Follow the detailed Directions for Use below. No other collection cassettes should be used with this assay. Oral fluid collected at any time of the day may be used.

[MATERIALS]

- | | |
|---|--|
| Materials Provided | Materials Required But Not Provided |
| <ul style="list-style-type: none"> Test Cassettes Security Seals Timer | <ul style="list-style-type: none"> Collectors Package Insert |
| <ul style="list-style-type: none"> Collection Chambers | |

[DIRECTIONS FOR USE]

Allow the test cassette, oral fluid collector, collection chamber to reach room temperature (15-30°C) prior to testing. Instruct the donor to not place anything in the mouth including food, drink, gum or tobacco products for at least 10 minutes prior to collection.


1. Sample Extraction

Take the oral fluid collector out of the sealed pouch and insert sponge end into the mouth. Actively swab the inside of the mouth for at least 3 minutes to collect oral fluid until the sponge becomes fully saturated. Gently press the sponge between the tongue and teeth to assist saturation. No hard spots should be felt on the sponge when saturated.

Note: Full saturation of sponge is very important to get enough oral fluid for processing. Otherwise, false or invalid results may be observed.

Take the collection chamber out of the sealed pouch and have it unscrewed. Remove the oral fluid saturated collector from the mouth and place it into the collection chamber. Press sponge fully against the strainer to release oral fluid to the arrowed point (level) and then screw it back tightly. Discard the collector.


Gently shake the collection chamber for 10 seconds and then leave it for 10 minutes.

2. Test Reaction

- Take test cassette out from sealed pouch, place it on a clean and level surface.
- Again gently shake the collection chamber for 10 seconds and snap the tip open, invert it and then transfer 3 drops of oral fluid specimen (approximately 120 µL) into specimen well(s) of the test cassette. Start the timer.
- Wait for the colored line(s) to appear. Read results at 5 minutes. Do not read results after 15 minutes.
- Apply a security seal over screw cap of the collection chamber and send to the laboratory for confirmation if necessary.

Step 3

[INTERPRETATION OF RESULTS]

(Please refer to the previous illustration)

NEGATIVE: Two colored lines appear. One colored line should be in the control region (C) and another colored line should be in the test region (T). This negative result indicates that the Marijuana concentration is below the detectable level of 3.5 ng/mL.

NOTE: The intensity of the color in the test line region (T) may vary depending on the concentration of Marijuana present in the specimen. Therefore, any shade of color in the test line region (T) should be considered negative.

POSITIVE: One colored line appears in the control region (C). No line appears in the test region (T). This positive result indicates that the Marijuana concentration is above the detectable level of 3.5 ng/mL.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test cassette immediately and contact your local distributor.

[QUALITY CONTROL]

A procedural control is included in the test. A colored line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

[LIMITATIONS]

- The THC Rapid Test Cassette (Oral fluid) provides only a qualitative, preliminary result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.
- A positive result indicates presence of the drug or its metabolites but does not indicate the concentration of drug in the specimen or the route of administration.
- A negative result may not necessarily indicate drug-free specimen. Negative results can be obtained when drug is present but below the cut-off level of the test.

[PERFORMANCE CHARACTERISTICS]

Accuracy

A side-by-side comparison was conducted using the THC Rapid Test Cassette and GC/MS at the cut-off of 3.5ng/mL. Testing was performed on 250 clinical specimens previously collected from subjects present for Drug Screen Testing. The following results were tabulated:

Method	GC/MS			Total Results
	Results	Positive	Negative	
THC Rapid Test Cassette	Positive	84	4	88
	Negative	2	160	162
	Total Results	86	164	250
	% Agreement	97.7%	97.6%	97.6%

Analytical Sensitivity

A Phosphate-buffered saline (PBS) pool was spiked with drugs to target concentrations of $\pm 25\%$ cut-off, $\pm 50\%$ cut-off and 300% cut-off and tested with the THC Rapid Test Cassette. The results are summarized below.

THC Concentration (ng/mL.)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0	30	30	0
1.75	-50%	30	30	0
2.625	-25%	30	27	3
3.5	Cut-off	30	12	18
4.375	+25%	30	8	22
5.25	+50%	30	0	30
10.5	300%	30	0	30

Analytical Specificity

The following table lists compounds and their respective concentrations in oral fluid that yield a positive result in the THC Rapid Test Cassette (Oral fluid) at 5 minutes.

Compound

Compound	Concentration (ng/mL)
11-nor- Δ^9 -THC-COOH	3.5
11-nor- Δ^9 -THC-COOH	2.5
Δ^9 -THC	40
Δ^9 -THC	20
Cannabidiol	2,000

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free oral fluid or Marijuana positive oral fluid. The following compounds show no cross-reactivity when tested with THC Rapid Test Cassette (Oral fluid) at a concentration of 100 µg/mL.

Non Cross-Reacting Compounds			
4-Acetamidophenol	Deoxycorticosterone	Dextromethorphan	(+)-3,4-Methylenedioxy-Prednisolone
Acetophenetidin	Diazepam	amphetamine	Prednisone
N-Acetylprocainamide	Diclofenac	methamphetamine	Promazine
Acetylsalicylic acid	Diltiazem	Methylphenidate	Promethazine
Aminopyrine	Digoxin	Methyprylon	D,L-Propanolol
Amtryptiline	Diphenhydramine	Morphine-3- β -D-glucuronide	D-Propoxyphene
Amobarbital	Doxylamine	Nalidixic acid	D-Pseudoephedrine
Amoxicillin	Ergonine	Nalorphine	Quinidine
	hydrochloride	Naloxone	Quinine
Ampicillin	Ergometrine	Naltrexone	Ranitidine
L-Ascorbic acid	(-)- α -Ephedrine	Naproxen	Salicylic acid
D,L-Amphetamine	Erythromycin	Nicardipine	Secobarbital
L-Amphetamine	β -Estradiol	Nisoldipine	Serotonin
Apomorphine	Aspartame	Oestrone-3-sulfate	(5-Hydroxy)amine
	Atropine	Estrone-3-sulfate	Sulfamethazine
Benzoic acid	Benzodiazepine	Ethynodiolide	Sulindac
Benzolic acid	Chloralose	Fenoprofen	Tetrazepam
Benzoyllecgonine	Chlorotetracycline	Furosemide	Tetracycline
Benzphetamine	Cinchonidine	Gentisic acid	Tetrahydrocortisone
Bilirubin	Cinchonine	Hemoglobin	Tetrahydrocortisone, 3-Acetate
(\pm)-Brompheniramine	Hydrochlorothiazide	Hydralazine	Tetrahydrocortisone
Caffeine	Hydrocodone	Isoxsuprime	Thiamine
Cannabidiol	Hydrocortisone	Ketamine	Thioridazine
Chlorhydrate	O-Hydroxyhippuric acid	Ketoprofen	D,L-Thyroxine
Chloramphenicol	3-Hydroxytyramine	Labetalol	Triamterene
Chlordiazepoxide	Ibuprofen	Levorphanol	Trifluoperazine
Chlorothiazide	Imipramine	Loperamide	Trimethoprim
(\pm) Chlorpheniramine	Iproniazid	Maprotiline	Trimipramine
Chlorquine	(\pm) - Isoproterenol	Meprobamate	Tryptamine
Cholesterol	Isoxsuprime	Methadone	D,L-Tryptophan
Clomipramine	Ketamine	Methoxyphenamine	Tyramine
Clonidine	Ketoprofen	Phenobarbital	D,L-Tyrosine
Cocaine	Labetalol	Phentermine	Uric acid
hydrochloride	Levorphanol	Phenylpropanolamine	Verapamil
Codine	Loperamide	Phenothiazine	Zomepirac
Cortisone	Maprotiline	Phenyltolamide	
(-)- Coline	Meprobamate	β-Phenylethylamine	
Creatinine	Methadone	Methoxyphenamine	
		Phenylpropanolamine	

[BIBLIOGRAPHY]

- Schramm, W. et al, "Drugs of Abuse in Oral fluid: A Review," J Anal Tox, 1992 Jan-Feb; 16 (1), pp 1-9
- Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488

Index of Symbols

	Consult instructions for use or consult electronic instructions for use		Contains sufficient for <n> tests		Temperature limit
	In vitro diagnostic medical device		Batch code		Catalogue number
	Authorized representative in the European Community/European Union		Use-by date		Do not re-use
	Do not use if package is damaged and consult instructions for use		Manufacturer		Caution

Hangzhou AllTest Biotech Co.,Ltd.

#550 Yizhe Street

Hangzhou Economic & Technological Development Area

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MedNet EC-REP GmbH